Claristep®

Syringeless, reliable, ultrapure filters from Sartorius

Claristep® Filters

Fill the liquid to be filtered into the liquid reservoirs of the Claristep® Filter, preferably with a pipette. The liquid filling volume should be at least 60 µl and not more than 600 µl. The recommended fill level is 400 µl (c).

Now close the lid of the Claristep® Station. This will close and press down the piston head of the Claristep® Filters (d).

Make sure to apply uniform pressure when closing the lid (e).

After closing the lid completely, hold it for three seconds and open it again. The Claristep® Filters are now fully closed and the filtrate is in the collection vials. Now remove the Claristep® Filters and dispose of them properly. The sample vials can now be removed and sealed if required (f).

Caution!

- For optimum handling of Claristep® Filters, we recommend using the Claristep® Station, including the Tray.
- If the Multiplex-Station is not used and the Claristep® Filter is manually closed, we recommend the use of a Claristep® Tray, since the supporting edge for the filters prevents collection vessels from fracturing (glass breakage).
- If neither the Multiplex-Station nor the Tray are used, there is no guarantee that implementation will be optimal.

Note

1. The membrane in Claristep® Filters is securely sealed from above to the edge of the Claristep® Filter's filter reservoir without any additional cover. This guarantees a very low dead volume of only approx. 30 µl, and small sample volumes of 60 µl can also be filtered. The visible sealing edge is not a quality limitation.


2. The piston head of the Claristep® Filter is equipped with a spacer that touches down on the sealed seam when the lid is closed and protects the membrane from mechanical damage.

3. One person should operate the Multiplex-Station by hand, with normal effort.

4. Use a Claristep® Filter for only one sample to reliably prevent carrying over residues from one sample to the next.

5. You will find an overview of chemical compatibilities in the Claristep® Chemical Compatibility Guide on our website.

Claristep® Station incl. replaceable Claristep® Tray

For optimal filtration with Claristep®, we recommend using the Claristep® Multiplex-Station incl. the replaceable Claristep® Tray (Order No.: 17C-M8).

Claristep® Single Tray

The Tray is also available separately (Order No.: 17C-S1).

The Claristep® Multiplex-Station and Single Tray are not supplied with Claristep® single-use filters and can be ordered separately using the specified article numbers.

Applications

Claristep® for particle removal from aqueous and organic liquids.

Instructions for Use

Place suitable collection vials (12 x 32 mm vials) in the Claristep® Tray (a).

Then place the Claristep® Filter above the collection vial. The leading edge of the Filter and supporting edge in the Tray ensure proper positioning. In addition, the Claristep® Tray has a custom-fit supporting edge, so that the Claristep® Filter does not rest directly on the collection vial. This minimizes the risk that glass vials will break (b).

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Table 1

Specifications for Claristep®

Housing material	Polypropylene (PP)	
Membrane	RC = Regenerated Cellulose	
Filtration Area	Effective membrane filtration diameter 9.7 mm Ø = 0.74 cm² effective filtration area	
Hold up volume	30 µl	
Application limits	Max. recommended operating pressure for 400 µl volume 2.5 bar / 36 psi	
Housing burst pressure	> 7 bar 102 psi	
Max. temperature	121°C, 30 min (autoclavable)	
Sterilization	Non-sterile Claristep® can be autoclaved or sterilized by ethylene oxide sterilization (EO)	
Claristep® Type	RC 0.2 µm	RC 0.45 µm
Non-sterile packs: 96 or 480 units	96: #17C07FT---96 480: #17C07FT---480	96: #17C06FT---96 480: #17C06FT---480
Bubble point	with water ≥ 3.0 bar 44 psi	with water > 5.0 bar 73 psi
Water penetration point	-	-
Sterile filtration capability* acc. to BCT	no	no
Main applications	universal hydrophilic filter with very low non-specific binding, compatible with many aggressive hydrophilic solvents	

* According to bacterial challenge test (BCT) with 10^7 Brevundimonas diminuta.
Non-sterile RC Minisart types are optimized for sample preparation and are not suitable for sterile filtration according to the BCT. All other non-sterile Minisart types listed above can be sterilized by autoclaving or EO before use for sterile filtration.

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
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Specifications subject to change without notice. Printed and copyrighted by Sartorius Stedim Biotech GmbH
Publication No.: SL-6207-p151001
Ver. 10 | 2015