

Tube Mill control Disposable system

IKA®



designed
for scientists

HUBERLAB. AG
Industriestrasse 123
4147 Aesch

T 061 717 99 77
F 061 711 93 42

www.huberlab.ch
info@huberlab.ch

HUBERLAB.

committed to science

Tube Mill control | Disposable System

IKA® introduces the world's first disposable grinding system for safe, instant and precise milling results. Its unique and compact design makes the unit space saving and ultra-portable. The disposable grinding chamber eliminates the possibility of cross-contamination and saves you cleaning costs and time.

The Tube Mill control is a batch mill for grinding soft, fibrous, hard and brittle materials (Mohs hardness up to 5). The transparent grinding chamber and cover facilitate observation at all times. Convenient and safe to use while assuring high safety and reproducibility to cover a broad range of applications. Amongst other applications, the mill is suitable for grinding seeds, such as corn and wheat. The ability to cool the sample with dry ice expands applications tremendously. During development of the mill, particular emphasis was placed on user safety.

The Tube Mill control is the world's first patented batch mill with disposable grinding chamber, designed and manufactured exclusively by IKA®.

IKA®+

World's first



- > Disposable grinding chamber
- > Cross-contamination free
- > No cleaning required
- > Precise milling results
- > Large application range



3 Year warranty*

* 2+1 years after registering at www.ika.com/register, glassware and wearing parts excluded



reddot design award
winner 2013



German
Design Award
NOMINEE 2014

Protection class according to DIN EN 60529: IP 30





Applications & Industries

The Tube Mill control is a highly-versatile milling device suitable for a broad range of applications used in various industries



> Food

- Rice
- Wheat
- Cobnut
- Coffee beans
- Spices
- Seeds
- Parsley
- Chocolate
(with dry ice)

> Pharmacy

- Vitamin tablets
- Tea leaves
- Pastilles
(with dry ice)
- Glauber salt
- Salt of hartshorn
- Blond plantain
- Sour orange paring
- Hawkbite roots
- Calamus roots



> Cosmetics

- Color pigments
- Rubber benzoe
- Bees wax
(with dry ice)



> Chemical Industry

- Rubber
- PE
- PET flakes
- Molecular sieve

> Biology

- Leaves (with dry ice)
- Grass (with dry ice)
- Tobacco
- Fresh ginger
(with dry ice)



> Medicine / Forensic

- Chicken bones
- Chicken skin
(with dry ice)
- Pig claw
(with dry ice)
- Chicken gristle
- Teeth
- Bone



> Renewable energy

- Straw
- Wood pellets
- Woodwool + wax
- Solid recovered fuel pellet
- Tetrapack
- Charcoal
- Chipped wood



> Building Materials Industry

- Gypsum
- Marble

Cryo grinding for advanced results

Typical applications for sample embrittlement with dry ice:

- > chocolate
- > bread
- > nuts
- > soil samples
- > gummy bears
- > leaves
- > meat
- > sausages
- > some plastics
- > beef
- > bones
- > feedstuffs
- > tobacco
- > grass

The Tube Mill control can also be used to process moist, fatty, elastic and fibrous samples.

Dry ice is introduced directly into the grinding chamber in order to embrittle the sample. The insulating effect of the plastic chamber allows minimal amounts of dry ice to be used. The cold remains in the milling chamber, allowing the user to handle the grinding chamber even after cooling. This greatly increases the range of applications for which the Tube Mill control can be used.



Tube Mill control | Technical data

**World's
first**



| Technical data | Tube Mill control |
|---|---------------------------------------|
| Process type | batch |
| Operating principle | cutting / impact |
| Motor rating input / output | 100 / 80 W |
| Speed range | 5000 – 25,000 rpm |
| Max. circumferential speed | 65 m/s |
| Max. usable volume | 40 ml |
| Timer | 5 s – 3 min |
| Interval timer | 5 – 60 s |
| Display | OLED |
| Max. Feed hardness | 5 Mohs (manganese or apatite: 5 Mohs) |
| Max. granularity of task | 10 mm |
| Mill feed can be cooled in milling chamber with dry ice | yes |
| Dimensions (W x D x H) | 180 x 300 x 170 mm |
| Weight | 2.7 kg |
| Permissible ambient temperature | 5 – 40 °C |
| Permissible relative moisture | 80% |
| Protection class according to DIN EN 60529 | IP 30 |
| USB interface | yes |
| Voltage | 220 – 240 V |
| Frequency | 50/60 Hz |

Ident. No. 0004180000*

* Two single grinding chambers are included in the delivery

IKA®+

Special safety features

- > The mill can only operate if the hood is closed
- > The motor does only operate with a correct grinding chamber
- > The system recognizes if the grinding chamber is not properly closed and the machine will not operate
- > The grinding chamber cannot be opened during the process
- > The motor is fitted with a labyrinth seal, preventing dust from entering the motor



Disposable grinding chamber, 40 ml

| | Ident. No. |
|---|------------|
| MT 40.10 | 0004425000 |
| Disposable grinding chamber 40 ml (10 pieces/pack) | |
| MT 40.100 | 0020001173 |
| Disposable grinding chamber 40 ml (100 pieces/pack) | |
| MTC 40.100 | 0020001182 |
| Cover for MT 40 | |
| MT 40.10 steril | 0020001984 |
| Disposable grinding chamber steril 40 ml (10 pieces) | |
| MT 40.100 steril | 0020001985 |
| Disposable grinding chamber steril 40 ml (100 pieces) | |

The disposable grinding chamber with a stainless steel beater reduces soft, medium, hard and brittle materials with a Mohs hardness of up to 5 (manganese or apatite: 5 Mohs). The chamber is made of transparent plastic so grinding tests can be observed at any time. In addition the chamber shows excellent resistance to chemicals and temperature.



Multiple grinding chamber, 40 ml

| | Ident. No. |
|--|------------|
| MMT 40.1 | 0020003165 |
| MMT 40.1 Stainless Steel package, includes 1 MMT 40 chamber, 25 sealings, 5 beaters, 5 couplings | |
| A-MMT 40.100 | 0020003378 |
| A-MMT 40.100 Abrasion Set Spare Parts, includes 100 sealings, 10 beaters, 10 couplings | |



Disposable grinding chamber:
 Stops cross-contamination
 No cleaning required
 Easy handling

100 ml chamber coming soon.
 Please contact your IKA® representative for further information

The Multi-use Milling Tube MMT 40.1 with a maximum volume of 40 ml can be used and cleaned in a dishwasher multiple times. The package includes a rich set of spare parts so that wearing parts can be replaced if necessary. Therefore it is possible to carry out a number of experiments with the grinding chamber depending on the nature of the sample.

Tube milling | Smooth process and easy storage



Step 1 | Fill the sample in the grinding chamber



Step 2 | Attach the grinding chamber onto the Tube Mill



Step 3 | Start the milling process



IKA® +

No cleaning costs



The grinding chamber can either be disposed of after the test or it can be used for storage of the processed sample. This new procedure will save on both time and money. As no cleaning of the tube system is required, the user is safe from aerosol formation that frequently occurs during cleaning procedures.

After grinding, a part of the sample will be analyzed. The remaining sample can either be discarded or it can be stored as a reference sample directly in the grinding chamber. In the later case, grinding chambers can be labeled and either stored in a refrigerator or in a drying room. Reference samples can be re-analyzed and traced at any time.



Step 4 | Grinding the sample



Step 5 | Remove the grinding chamber



Step 6 | Remove the grinding sample



IKA® offers more



Worldwide service network – direct contact in your region

Our dedicated team of engineers provides comprehensive worldwide technical service. Please feel free to contact IKA® directly or your dealer in case of any service questions.

For spare parts IKA® guarantees 10 years of availability. In the event of an equipment malfunction or technical questions regarding devices, maintenance and spare parts, please call us at 00 8000 4524357 (00 8000 IKAHELP) or send an email to service@ika.com



IKA® Application Support

Our Application Center spans 400 sqm and offers modern facilities for presenting and testing lab devices and processes. This brings us even closer to our customers and improves our service. Here, prospective buyers and customers can test processes that involve stirring, shaking, dispersing, grinding, heating, analyzing and distilling.

Call us at 00 8000 4522777 (00 8000 IKAAPPS) or send an email to applicationsupport@ika.com or visit our website at www.ika.com/applicationsupport



Customizing Center

It is important that IKA® products work for your application. We have a special program: product solutions tailored to your needs.

Should you not find the appropriate device in our standard product range, please send us your requested specifications through the online form. Our team will determine its feasibility and offer a solution to you.

Please visit www.ika.com/customizingcenter to review already implemented product modifications.

Tube Mill control | Your benefits




Adjustable safety speed and time




Interval operation available



USB interface to control and document all the parameters and for updating your firmware



Special safety features



Clearly arranged, multi-lingual OLED display



Quiet operation



3 Year warranty*

* 2+1 years after registering at www.ika.com/register, glassware and wearing parts excluded

Protection class according to DIN EN 60529: IP 30



Please visit
www.ika.com

Service | FAQ

Is it possible to use the grinding chamber more than one time?

We recommend to use the grinding chamber only once to avoid cross-contamination

What about the cleaning methods of the grinding chamber?

Before the first use, the grinding chamber can be autoclaved

What material are the grinding chamber, knife and vials made of?

The grinding chamber is made of PP, the knife is made of spring steel 1.4310 and vials are made of PA

Can standard grinding chamber be used with dry ice?

Yes, the grinding chamber can be used with dry ice

What about the end fineness of samples?

The end fineness is between 1 – 100 μm
(depends on sample)

What about the minimum quantity for the grinding chamber?

One corn

What about the Mohs hardness of the samples?

The maximum Mohs hardness for samples is 5

