



# neoform<sup>®</sup> Rapid



## Rapid disinfectant for medical devices and surfaces

### Ready-for-use solution

### Fields of application:

Specific and rapid disinfection of

- surfaces of medical devices, e.g. surgical tables and other surfaces in hospitals and doctors' and dentists' practices
- surfaces and equipment in the pharmaceutical and cosmetics area
- surfaces in the laboratory area, e.g. laboratory benches

Also suitable for use with the Dr. Weigert fleece wipe dispensing system neoform wipes RTF.

### Performance spectrum:

- Bactericidal, yeasticidal, fungicidal and tuberculocidal activity confirmed in accordance with VAH<sup>1</sup> and EN methods<sup>2</sup>; activity against enveloped viruses confirmed according to the test guidelines of the RKI/DVV<sup>3</sup> by certification
- Also limited spectrum virucidal activity (active against noro- and adenovirus) as well as active against rotavirus
- The application in the Dr. Weigert fleece wipe dispensing system neoform wipes RTF fulfils the requirements of the current VAH recommendation<sup>4</sup>
- The surfaces to be disinfected must be compatible with alcoholic disinfectants
- Not suitable for disinfecting acrylic glass (perspex)
- VAH<sup>1</sup>-listed
- Included in the IHO<sup>5</sup> disinfectant list

### Special properties:

- Wide activity spectrum, even with short contact times
- Rapid and stain-free drying, no rinsing with fresh water necessary
- Free of aldehydes, colourants and perfumes

### Application and dosage:

Apply neoform Rapid undiluted on clean and dry surfaces and disinfect by wiping. Surfaces must be completely wetted with neoform Rapid. To achieve the desired activity, the below-mentioned contact times must be strictly adhered to.

For instructions on how to use the product with the neoform wipes RTF fleece wipe dispensing system please refer to the neoform wipes RTF data sheet. The soaked neoform wipes can be used for up to 28 days when the wipes system is kept tightly closed.

	Clean conditions, 20 °C (cleaned surfaces)
bactericidal (VAH <sup>1</sup> , EN 13727, EN 16615)	undiluted, 1 min
yeasticidal (VAH <sup>1</sup> , EN 13624, EN 16615)	undiluted, 1 min
fungicidal (EN 13697)	undiluted, 3 min
tuberculocidal (EN 14348)	undiluted, 30 sec
active against enveloped viruses (incl. HIV, HBV, HCV) (EN 14476, RKI/DVV <sup>3</sup> )	undiluted, 15 sec
limited virucidal activity (adenovirus, norovirus, incl. all enveloped viruses) (EN 14476)	undiluted, 2 min
active against norovirus (MNV <sup>6</sup> ) (EN 14476)	undiluted, 2 min
active against adenovirus (EN 14476)	undiluted, 1 min
active against rotavirus (EN 14476)	undiluted, 1 min
VAH <sup>1</sup> -list	undiluted, 5 min

Due to its alcohol content, neoform Rapid must not be applied in too large amounts (max. 50 ml per m<sup>2</sup> on the surfaces to be treated and max. 100 ml per m<sup>2</sup> base area). Please thoroughly read and observe the hazard and precautionary statements!

### Notes on application:

- For professional use only.



# neoform<sup>®</sup> Rapid

- For economical and controlled dosage the use of manual dosing aids is recommended. Please contact us.
- It is generally recommended to wear gloves when working with disinfectants.
- Reprocessing should comply with all ordinances pursuant to the Medical Device Directive and should be performed with appropriate validated processes
- Please observe the reprocessing recommendations given by the manufacturer of the medical devices and surfaces.
- Do not mix with other products.
- *Use disinfectant safely. Always read the label and product information before use.*

## Expert reports:

The disinfecting activity has been confirmed by certification. We will be pleased to provide certificates on request.

## Technical data:

Viscosity	< 10 mPa s (undiluted, 20 °C)
Density	approx. 0.9 g/cm <sup>3</sup> (20 °C)
Flash point	approx. 27 °C according to DIN 22719


## Ingredients:

Active substances in 100 g:  
35.0 g propan-1-ol  
25.0 g ethanol

CE-mark:  0297

neoform Rapid complies with European guidelines for medical devices.

## Storage information:

Always store at a temperature between -15 °C and 25 °C. Keep container tightly closed. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see safety data sheets. These are available at [www.drweigert.com](http://www.drweigert.com)

under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet

MB 5643/3-2 Date of Issue: 07/2019

- 1 Verbund für Angewandte Hygiene [Association of Applied Hygiene, Germany]
- 2 European test methods
- 3 according to the guidelines of the Robert Koch-Institute (RKI) and the Deutschen Vereinigung zur Bekämpfung von Viruskrankheiten (DVV) [German Association for the Control of Virus Diseases]
- 4 Empfehlung zur Kontrolle kritischer Punkte bei der Anwendung von Tuchspendersystemen im Vortränkssystem für die Flächendesinfektion (Hyg Med 2012; 37-11) [Recommendation for checking critical points when using wipes dispenser systems in a pre-soak system for surface disinfection]
- 5 Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]
- 6 Murine norovirus

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.

**HUBERLAB.**

committed to science

HUBERLAB. AG  
Industriestrasse 123  
4147 Aesch

T +41 61 717 99 77  
F +41 61 711 93 42

[www.huberlab.ch](http://www.huberlab.ch)  
[info@huberlab.ch](mailto:info@huberlab.ch)