

# Cryo.s™ Rack Scanner

Ideal for Labs Scanning Frozen Racks

## Key Facts

- ☞ CCD based scanner for 2D coded tubes and SBS-formatted racks
- ☞ Scans 1 rack in 8-10 sec
- ☞ Active and passive condensation reduction
- ☞ Single tube scan function
- ☞ One touch scan button
- ☞ Robust metal housing
- ☞ Single to use scanning software with customizable output
- ☞ Toolkits (SDK, ActiveX) for integration into sample management software\* or LIMS\* as part of software package



## Ordering Information

Cat. No.	Product Description	Box contents
849 050	Cryo.s™ Rack Scanner, US edition	1 scanner, power cables, USB cables, drivers and manual
849 070	Cryo.s™ Rack Scanner, EU/UK/Asia edition	1 scanner, power cables, USB cables, drivers and manual

\*A sample management software is not part of the included software package.