

# New Plate Printing

Article No. 657160

## CELLSTAR® Multiwell Plate, 6 Well, PS, clear, TC

In addition to the established design and proven quality the above mentioned plate will come in future with a black printing as novel feature.

The printing was designed to improve the identification of the respective plate as well as its usability. It includes:

- **Alphanumeric coding for faster orientation**
- **Type of surface treatment for simple and fast plate identification**
- **Lot number and expiration date**  
(For technical reasons, bar coded plates can just be printed with alphanumeric coding.)

Furthermore, a new frosted field at the side of the plate enables individual labelling by the customer.

The existing article number **657160** will be retained for the new plates. For a smooth transition, boxes containing the printed version are labeled with **“New Design”**.



6 well multiwell plate with black printing of:

- Alphanumeric coding
- Lot number
- Expiry date
- Surface treatment