

Ziath Handheld 2D Barcode Scanner



Features

- Intelligent and fully portable 2D barcode reader for single tubes
- Full day operation without the need to re-charge the battery
- 4 operating modes
- Cryoprotection provided as standard
- Up- or down-load data from host computer

Four Operating Modes

<p>1) Simple barcode scanning Display the barcode value from a scanned tube</p>	<p>3) Pick List Upload a list of samples to the Handheld and check them off one at a time</p>
<p>2) Check sample data Display complete sample details from a pre-loaded check list</p>	<p>4) Add sample data Edit an existing or add a new sample</p>



- **Dimensions = 155mm x 65mm x 38mm (length x width x max depth)**
- **Weight = 150g**
- **Battery charger input = 100-240VAC 60/50Hz, 260mA**
- **Battery charger output = 5VDC, 2.4A**
- **Battery life = Typically 8 hours between charges**
- **Data storage on instrument = Typically up to 60 000 samples**
- **Display size = 55mm x 90mm**



Included in the Handheld Package

- **Handheld instrument**
- **Battery charger**
- **Software package to provide user friendly interface between computer and instrument**
- **USB lead to connect instrument to host computer**

Ordering Information

Code	Description	Quantity
ZTS-HND	Ziath Intelligent Handheld 2D barcode reader	1