

# Model LD 79 Test-tube-rotator



- Provides mixing for diffusion, dialysing, extraction etc.
- Tube diameter range 4 - 40 mm
- Accepts all disks mentioned in the listing on front page
- Mixing action can be setted by:
  - speed control 5-50 rpm
  - angle of rotation 90° to 27°
- Quick/easy exchange of disks by unscrewing central nut
- Dual or triple disk support may be carried out
- Tubes may be exchanged while the rotator is running
- Dimensions:  
200(l) x 200(w) x 250(h) mm
- Weight 3,5 kg

50/60 Hz Cat.No.: 79000 (230V) / 79010 (115V)

## Electronically controlled stepper-motor provides:

- Digital read-out for - speed 5-50 rpm  
- time 0-99 hrs/1-59 min
- Electronic timer for "auto-stop"
- "Memory-function" for repeated operation.
- "Re-start-function" for continued operation after line-voltage failure including auto-stop
- "Settings saving" available

