

# Technical Datasheet

Date: 21.03.18

Page 1 of 2

EN ISO certified		
<b>Article - No:</b>	<b>11-702-1010</b>	<b>Revision: 1</b>
<b>Remarks</b>	<b>suitable for formalin and chloroform</b>	<b>single use only</b>

## Specification

<b>Product</b>	Screw cap container PP, 500ml, Ø90x91 mm, with fitted screw cap PP, with PE-pad, white, autoclavable up to 121°C (without pad)
<b>Dimensions</b>	container:
	diameter: 89,20 mm (+/- 0,2 mm)
	height: 90,50 mm (+/- 0,2 mm)
	cap:
	diameter: 96,00 mm (+/- 0,2 mm)
	height: 15,70 mm (+/- 0,2 mm)
	total height container and cap: 91,00 mm (+/- 0,2 mm)
<b>Colour</b>	white
<b>Nominal Capacity</b>	500ml
<b>Material Characteristics</b>	container: polypropylene (PP)
	cap: polypropylene (PP)
	pad: polyethylene (PE)-foam (400 kg/cbm)
<b>Temperature range</b>	-80°C - 70°C
<b>Autoclavability</b>	autocl. up to 121°C (without pad)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	flexible and soft
<b>Shelf - Life</b>	not limited at controlled storage conditions
<b>Packaging</b>	pcs per cs: 72
	pcs per carton: 72
	pcs per pallet: 1.728
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags&cases orig. closed