

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	02-401-0000	Revision: 2
Remarks	Fitting plug: 02-633-0000	single use only

Specification

Product	Test tube PS, conical bottom, 4ml, Ø11x75 mm, transparent, max. RCF 3.000g
Dimensions	outer diameter: 11,00 mm (+/- 0,2 mm)
	height: 75,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	4ml
Material Characteristics	polystyrene (PS)
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 100
	pcs per carton: 4.000
	pcs per pallet: 80.000
Minimum order qty.	4.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	