

Technical Datasheet

Date: 17.12.18

Page 1 of 1

EN ISO certified		
Article - No:	02-412-0001	Revision: 3
Remarks	Fitting plug: 02-633-0000	single use only

Specification

Product	Test tube PP, conical bottom, 4,5ml, Ø12x75mm, transparent, max. RCF 5.500g, autocl. up to 121°C
Dimensions	outer diameter: 12,00 mm (+/- 0,2 mm)
	height: 75,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	4,5ml
Material Characteristics	polypropylene (PP)
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 5.500g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 4.000
	pcs per pallet: 96.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department