

# Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>02-421-0000</b>	<b>Revision: 1</b>
<b>Remarks</b>	<b>Fitting plug: 02-673-0000</b>	<b>single use only</b>

## Specification

<b>Product</b>	Test tube PS, conical bottom, 12ml, Ø16x105 mm, transparent, max. RCF 3.000g
<b>Dimensions</b>	outer diameter: 16,00 mm (+/- 0,2 mm)
	height: 105,00 mm (+/- 0,2 mm)
<b>Colour</b>	transparent
<b>Nominal Capacity</b>	12ml
<b>Material Characteristics</b>	polystyrene (PS)
<b>Centrifugation</b>	max. RCF 3.000g
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
<b>Packaging</b>	pcs per bag: 1.000
	pcs per carton: 2.000
	pcs per pallet: 40.000
<b>Minimum order qty.</b>	2.000 pcs
<b>Storage</b>	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked

Approved by nerbe plus quality department