

Technical Datasheet

Date: 17.02.17

Page 1 of 1

		EN ISO certified
Article - No:	02-531-6007	Revision: 3
Remarks		single use only

Specification

Product	Test tube PS with fitted screw cap PE, 25ml, Ø25x90 mm, conical bottom, skirted, tube transparent, cap white, sterile R, max. RCF 4.000g, CE/IVD
Classification	CE/IVD
Dimensions	tube: outer diameter: 25,00 mm (+/- 0,2 mm) height: 90,00 mm (+/- 0,2 mm) cap: outer diameter: 31,00 mm (+/- 0,2 mm) height: 17,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap white
Nominal Capacity	25ml
Material Characteristics	tube: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-10°C – 70°C
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 100 pcs per carton: 500 pcs per pallet: 12.000
Minimum order qty.	500 pcs
Storage	dry condition: 40%-60% rel. humidity bags & cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	