

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	02-552-1000	Revision: 1
Remarks		single use only

Specification

Product	Test tube PP with enclosed screw cap PE, 30ml, Ø25x107 mm, conical bottom, skirted, tube transparent, cap natural, graduated, max. RCF 6.000g
Dimensions	outer diameter: 25,00 mm (+/- 0,2 mm) height: 107,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Graduation	graduated
Nominal Capacity	30ml
Material Characteristics	tube: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-20°C - 60°C
Centrifugation	max. RCF 6.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 500 pcs per carton: 500 pcs per pallet: 12.000
Minimum order qty.	500 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	