

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	08-022-2004	Revision: 1
Remarks		single use only

Specification

Product	Universal specimen container with spoon PP, 15ml, Ø20x76 mm, with fitted screw cap PE, skirted, container transparent, cap brown, graduated, max. RCF 2.500g, CE/IVD
Classification	CE/IVD
Dimensions	diameter: 20,00 mm (+/- 0,2 mm) height (without screw cap): 76,00 mm (+/- 0,2 mm)
Colour	container transparent, cap brown
Graduation	graduated
Nominal Capacity	15ml
Material Characteristics	tube: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-20°C - 65°C
Centrifugation	max. RCF 2.500g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Packaging	pcs per bag: 1.000 pcs per carton: 1.000 pcs per pallet: 24.000
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	