

Technical Datasheet

Date: 17.05.17

Page 1 of 1

EN ISO certified		
Article - No:	04-072-0390	Revision: 3
Remarks		single use only

Specification

Product	Reservoir PP, 140ml, high profile, 12 channels, transparent, np pcr ready
Dimensions	ANSI/SLAS 1
Colour	transparent
Nominal Capacity	140ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
	free of endotoxins
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 5
	pcs per carton: 25
	pcs per pallet: 1.500
Storage	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department