

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-023-0100	Revision: 1
Remarks		single use only

Specification

Product	Screw cap HDPE, for mailing container 30x85 mm and 31x126 mm, natural
Colour	natural
Material Characteristics	high density polyethylene (HDPE)
Temperature range	-50°C – 70°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Packaging	pcs per bag: 500 pcs per carton: 500 pcs per pallet: 54.000
Minimum order qty.	500 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	