

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-062-1000	Revision: 1
Remarks		single use only

Specification

Product	Mailing container PE, 34x34x85 mm, with absorbing inlay, with enclosed cap PE, natural
Dimensions	length: 34,00 mm (+/- 0,2 mm)
	width: 34,00 mm (+/- 0,2 mm)
	height: 85,00 mm (+/- 0,2 mm)
Colour	natural
Material Characteristics	polyethylene (PE)
Temperature range	-50°C – 70°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 500
	pcs per carton: 500
	pcs per pallet: 10.000
Minimum order qty.	500 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked

Approved by nerbe plus quality department