

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-222-0100	Revision: 1
Remarks		single use only

Specification

Product	Mailing container PP, 50ml, Ø44x114 mm, with absorbing inlay, screw thread and print according to EN 829, without cap, natural, autoclavable up to 121°C
Dimensions	diameter: 44,00 mm (+/- 0,2 mm)
	height: 114,00 mm (+/- 0,2 mm)
Colour	natural
Nominal Capacity	50ml
Material Characteristics	container: polypropylene (PP)
	inlay: cellulose
Temperature range	-196°C - 121°C
Autoclavability	autoclavable up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 250
	pcs per carton: 250
	pcs per pallet: 5.500
Minimum order qty.	250 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	