

Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-871-0000	Revision: 1
Remarks		single use only

Specification

Product	Drosophila vial PS, 20ml, Ø20x60 mm, flat bottom, highly transparent
Dimensions	inner diameter: 20,00 mm (+/- 0,2 mm)
	outer diameter: 22,00 mm (+/- 0,2 mm)
	height: 60,00 mm (+/- 0,2 mm)
Colour	highly transparent
Nominal Capacity	20ml
Material Characteristics	polystyrene (PS)
Temperature range	-10°C – 70°C
Product Characteristics	biologically inert
	physiologically harmless
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 2.000
	pcs per carton: 2.000
	pcs per pallet: 40.000
Minimum order qty.	2.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	