

Product Specification

SELECT PLUS ECO

Disposable Latex Gloves

- Powder free
- White
- Double chlorinated
- Non-sterile



SAFETY MEETS QUALITY

Material

Natural High Grade Rubber Latex

Performance requirements

Characteristics	Inspection Level	Acceptable Quality Level
Dimensions	S-2	4.0
Physical Properties	S-2	4.0
Freedom from Holes	G-1	1.5
Major Visual defects	G-1	2.5
Minor Visual defects	G-1	4.0

Dimensions

	1401	1402	1403	1404
Number	1401	1402	1403	1404
Size	XS = 5-6	S = 6-7	M = 7-8	L = 8-9
Width	75 ± 5 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Thickness

Measuring Points	Single Wall (mm)
Cuff (distance 25 mm to the beaded edge)	0.09 ± 0.02
Palm (in the middle of the palm)	0.12 ± 0.02
Fingertip (distance 15 mm to top of middle finger)	0.13 ± 0.03

Physical Properties

Standard	Before Aging	After Aging
Tensile Strength (MPa)	min. 6.0	min. 6.0
Elongation at Break (%)	min. 650	min. 500
Median force at break (MPa)	min. 18.0	min. 14.0

Powder

No powder lubricant added

Protein content

< 50 µg/dm²

Cuff

Beaded cuff

Design

Ambidextrous

Surface

Micro-roughened palm

Color

White

Quality assurance

Comply the requirements of ASTM D3578, EN 455 part 1-4, Manufactured under ISO 9001:2008, ISO 13485:2012, according GMP

Packing Style

150 pcs per dispenser, 10 dispensers per carton

Hint

The above information apply to the standard version. Subject to change.