



DR. WEIGERT
Systematic Hygiene

Medicine | Home Care | Laboratories | Kitchens



neoform®
wipes RTF: the
Effective Wipes
Dispenser System
**with Impressive
Contents.**

Efficient surface disinfection for medicine,
home care, laboratories and kitchens.

The **All-round System** for Every Challenge.

Benefit from our expertise!

For over 100 years, Dr. Weigert has been setting standards for professional cleaning and disinfection – in hospitals, geriatric and care institutions, professional kitchens and laboratories. Dr. Weigert offers hygiene expertise from a single source for its customers in all sectors.

The independent family-run company with its own Research & Development department including laboratories and production facilities has an international reach.

In all its activities, Dr. Weigert always focuses on quality, customer proximity, efficiency and sustainability. With its customer-specific solutions, Dr. Weigert enjoys an outstanding reputation.

The right solution every time – in a practical package.

neoform® wipes RTF – the efficient disposable system for safe surface disinfection.

Disinfection and cleaning of small and large surfaces in medical, nursing-care, laboratory and kitchen areas are a key part of the daily hygiene routine.

And this is precisely where the danger lies – as it has been shown that mistakes are more likely to be made during routine activities. Inadequate surface disinfection has serious consequences: germs can spread because of their ability to survive, thus increasing the risk of cross-contamination and infections.

The innovative disposable wipes dispenser system **neoform® wipes RTF** from Dr. Weigert is a solution with established efficiency here.

The all-round system, which can be filled with the proven Dr. Weigert product solutions for surface disinfection, provides maximum safety in all areas: the medicine, home home, laboratory and kitchen area.

Practical handling combined with the guaranteed active-ingredient release and disinfecting action of the soaked wipes makes routine jobs easier for staff and plays a major role in protecting personnel and patients. This system combines Dr. Weigert's hygiene expertise with a practical solution for your staff's day-to-day work.



**neoform®
wipes RTF –**
a System Concept
That Also
**Impresses
in Field Tests.**



Infection Risks Successfully **Minimised** in Hospitals.



Protection for patients and medical staff.

neoform® wipes RTF – hygiene safety even in areas with a particularly high infection risk.

The **neoform® wipes RTF** disposable wipes dispenser system is particularly suitable for use in hospitals and in high-risk areas.¹ In conjunction with the corresponding Dr. Weigert disinfection solutions, the risk of cross-contamination from biofilm is significantly minimised, for example.

Studies demonstrate the benefits of ready-to-use wipes dispenser systems with the possibility of enhanced compliance²:

- Staff more prepared for disinfection.
- Activities performed more quickly.
- Minimised risk of germ transmission through disposable wipes.

¹ "The Disinfection Commission of the German Association of Applied Hygiene (VAH) advises that wipes dispenser systems that have to be reprocessed – with the exception of alcoholic agents – should not be used in areas with particular infection risks." Source: HygMed 2014

² "The value of ready-to-use disinfectant wipes: Compliance, employee time and costs – Source: American Journal of Infection Control, 2014."



Outstanding Support in **Hygiene** Matters.

Surface disinfection time reduced by 23%! A great help for nursing staff.

In home care, hygiene should not be compromised due to time pressure or a shortage of specialist staff.

Surface disinfection is a core part of the daily hygiene routine in home care in particular. As well as being clean and free of dust and contamination, the surroundings must be clear of any infection risk to residents and staff with regard to the bacteria count and the killing of pathogenic or potentially pathogenic micro-organisms. The innovative wipes dispenser system **neoform® wipes RTF** makes surface disinfection much more efficient – the time required to clean a patient's room is reduced by 23%. This time saving is also reflected in staff costs.³

The cost-benefit analysis for use of a disposable wipes system versus multiple-use wipes for surface disinfection presents a clear result: the study clearly shows that cost optimisation and a high standard of quality are entirely compatible with each other.⁴

³ "Cost-benefit analysis for the use of disposable versus multiple-use wipes for surface disinfection/Klinikum Nürnberg."
Source: HygMed 2008

⁴ "The value of ready-to-use disinfectant wipes: Compliance, employee time and costs – Source: American Journal of Infection Control, 2014."



With **neoform®**
wipes RTF, You Ensure
Clean Results
in the Laboratory.



Maximum hygiene is the basis for every test set-up.

Surface disinfection in the laboratory forms the basis for perfect, reliable analysis results.

The new disposable wipes system **neoform® wipes RTF** has also been proven to be a perfect system for surface disinfection in laboratory operations. Particularly in a laboratory environment, where analyses are performed on a daily basis, contamination by germs from prior analyses is possible.

For this reason, recurring and efficient routine disinfection is essential to the attainment of clean results. In conjunction with the **neoform® wipes RTF** system, as well as using the proven Dr. Weigert hygiene products for the laboratory, use of a disposable system also gives you the certainty that no germs are transmitted due to inadequately reprocessed wipes dispensers.



Hygienically
Clean Surfaces
**are a Recipe
for Success in
the Kitchen.**

Impress with great-tasting food and dazzling hygiene.

In food processing and commercial kitchens, it is particularly important to eliminate contamination and the transfer of viruses and bacteria.

Heat, moisture and any number of organic materials – what sounds like optimum requirements for multiplication of micro-organisms – are everyday conditions in professional kitchens. No wonder special attention must be paid to hygiene here.

Ultimately, a commercial kitchen is also the ideal place for inadvertently passing on germs to lots of people. Dr. Weigert has already set quality standards in this field with its hygiene products. With the new **neoform® wipes RTF**, we now offer you a system solution that is perfectly suited to use in the kitchen.



A totally convincing system:



Safe disposable system:

- The wipes are ready to use after a 30-minute contact time.
- Integrated fleece roll for a reduced contamination risk.
- Also suitable for high-risk areas such as neonatology and haematological oncology as recommended by the VAH.¹
- Eliminates error-susceptible reprocessing.
- Tamper-proof cap.



Reduced time & increased quality:

- Guaranteed active-ingredient release: proven disinfecting action of the soaked wipes for the maximum service life of 28 days under EN 16615.
- Saves time-consuming and cost-intensive reprocessing.



The all-rounder from Dr. Weigert:

- Wide range of applications.
- Stackable and space-saving.
- Tested compatibility with neoform® products as per EN 16615.
- Meets the requirements of the VAH recommendation on the use of wipes dispenser systems in combination with Dr. Weigert surface disinfectants.²
- Recyclable plastic.



RTF:
ready to fill

neoform® wipes RTF
8 wipe dispensers per l.
4213 95



Resealable container:

- High level of hygiene safety.
- Prevents drying out.
- Low contamination risk.



Dr. Weigert wipes:

- Easy to remove.
- Lint-free/suitable for operating theatres and laboratories.
- 115 high-quality wipes per roll.
- Highly absorbent.
- Outstanding tensile strength.



Optimally designed dispenser opening:

- Special teardrop shape makes tearing off and closing easier.



Stable & integrated handle:

- Ergonomically designed handle enables easy transportation even when full.

The best solution for all areas:



neoform® Rapid



neoform® Classic



neoform® K sprint



neoform® K plus

The following pages provide a detailed overview of the **neoform® products** that can be combined with **neoform® wipes RTF**.


neoform® Rapid & neoform® Classic for medicine, home care and laboratories.

neoform® Rapid¹

Fast-acting disinfectant for:

- Surfaces of medical devices
- Surfaces in hospitals, surgeries, dental practices and care institutions
- Surfaces and equipment in the pharmaceutical, cosmetics and laboratory sectors
- Bactericidal, yeasticidal, fungicidal, tuberculocidal action, active against enveloped viruses, limited spectrum of virucidal activity (effective against noroviruses² and adenoviruses) and effective against rotaviruses; confirmed by experts
- Included in the IHO³ list of disinfectants



12x 
neoform® Rapid
12 x 0.75 l bottle
564347

neoform® Classic¹

Amine-based disinfecting cleaner for:

- Surfaces of medical devices and medical equipment as well as other surfaces in hospitals and surgeries
- Surfaces in geriatric and care institutions, laboratory areas and public areas
- Bactericidal, yeasticidal and tuberculocidal action, active against enveloped viruses, limited spectrum virucidal activity (effective against noroviruses² and adenoviruses); confirmed by experts
- Included in the IHO³ list of disinfectants



 5 L
neoform® Classic:
5 l canisters
447035

6x  2 L
neoform® Classic:
6 x 2 l bottle
447042

Application:

Depending on the product, recommended for disinfecting cleaning or fast-acting disinfection of medical devices and surfaces in hospitals, surgeries, laboratories and care institutions as well as in the pharmaceutical and cosmetics sectors.

Suitable for application with **neoform® wipes RTF**. For correct preparation of the soaked wipes, please comply with the application recommendation. For safe disinfection, thorough wetting with the appropriate contact time must be ensured. Detailed application parameters can be found in the product data sheet (www.drweigert.de).

Application recommendation with neoform® wipes RTF:

neoform® Rapid

Bottle: 0.75 l

Wipes dispenser: 2.25 l (3 x 0.75 l)

Undiluted

neoform® Classic


Bottle: 2 l

Wipes dispenser: 2.5 l

Application solution
At least 0.5%

Our expertise –

- Compatibility/disinfecting action of the soaked wipes confirmed for the maximum service life of 28 days under EN 16615
- Stable wipes dispenser system allows practical handling and filling
- Disposable system



neoform® K sprint & neoform® K plus for professional kitchens and the food industry.

neoform® K sprint¹

Fast-acting disinfectant for:

- Surfaces in contact with food or products and other surfaces in the food industry and in professional kitchens, and also usable for disinfection in cooling areas
- Residue-free drying without rinsing, no colourants or perfumes
- Bactericidal, yeasticidal, fungicidal, tuberculocidal action, active against enveloped viruses, limited spectrum virucidal activity (effective against noroviruses² and adenoviruses) and effective against rotaviruses; confirmed by experts
- Included in the IHO³ list of disinfectants



neoform® K sprint:
10 l canisters
5644 30



12x
neoform® K sprint:
12 x 0.75 l bottle
5644 47

neoform® K plus¹

Disinfecting cleaner with outstanding cleaning performance for:

- Surfaces and equipment in the food industry and in professional kitchens
- Excellent material compatibility, good rinsability, no colourants or perfumes
- Bactericidal, yeasticidal and limited virucidal effectiveness, confirmed by experts
- Included in the IHO³ list of disinfectants



neoform® K plus:
200 kg drum
3223 13



neoform® K plus:
10 l canisters
3223 30



6x
neoform® K plus:
6 x 2 l bottle
3223 42

Application:

Depending on the product, for disinfecting cleaning or fast-acting disinfection of surfaces and equipment in the food industry and in professional kitchens as well as surfaces and equipment that do not come into contact with food, such as in the pharmaceutical and cosmetics industries.

Suitable for application with **neoform® wipes RTF**. For correct preparation of the soaked wipes, please comply with the application recommendation. For safe disinfection, thorough wetting with the appropriate contact time must be ensured. Detailed application parameters can be found in the product data sheet (www.drweigert.de).

Application recommendation with neoform® wipes RTF:

neoform® K sprint

Bottle: 0.75 l

Wipes dispenser: 2.25 l (3 x 0.75 l)

Undiluted

neoform® K plus

Bottle: 2 l

Wipes dispenser: 2.5 l

Application solution
At least 0.5%

Your benefits:

- Safe surface disinfection for patient and staff protection
- Possibility of enhanced compliance
- Elimination of time-consuming and cost-intensive reprocessing as well as the need to make up a stock solution every working day

¹ Use disinfectants safely. Always read the label and product information before use.

² Murine norovirus

³ Industrial Association for Surface Protection and Hygiene, Germany