

## Product Specification

# SOFT NITRIL BLUE PREMIUM

## Disposable Nitrile Examination Gloves

- Powder free
- Blue
- Latex free
- Double chlorinated
- Chemical resistance protective gloves
- Suitable for handling food
- Non-sterile



CE 0121



SAFETY MEETS QUALITY

### Material

100% Synthetic Nitrile Latex

### Performance requirements

Characteristics	Inspection Level	Acceptable Quality Level
Dimensions	S-2	4.0
Physical Properties	S-2	4.0
Freedom from Holes	G-1	1.5
Major Visual defects	G-1	2.5
Minor Visual defects	G-1	4.0

### Dimensions

Number	8020-374	8030-374	8040-374	8050-374
Size	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	≥ 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

### Thickness

Measuring Points	Single Wall (mm)
Cuff (distance 25 mm to the beaded edge)	0.08 ± 0.03
Palm (in the middle of the palm)	0.11 ± 0.03
Fingertip (distance 15 mm to top of middle finger)	0.13 ± 0.03

### Physical Properties

Standard	Before Aging	After Aging
Tensile Strength (N)	min. 6.0	min. 6.0
Elongation at Break (%)	min. 500	min. 400
Modulus of elasticity (MPa)	min. 14.0	min. 14.0

### Powder

No powder lubricant added

### Protein content

No proteins contained

### Cuff

Beaded cuff

### Design

Ambidextrous

### Surface

Micro-roughened palm

### Color

Blue

### Quality assurance

EG construction sample test IFA – CE 0121

### Performance standards

PSA – EN 374 CAT III, EN 388, EN 420  
 MPG – EN 455 part 1-4  
 EG-Standard 89/686/EWG  
 EG-Standard 93/42/EWG, Appendix II.3 + V.3  
 EN ISO 9001:2008  
 EN ISO 13485:2012 + AC:2007

### Packing Style

100 pcs per dispenser, 10 dispensers per carton

### Hint

The above information apply to the standard version.  
 Subject to change.