

# Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>02-841-2000</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Cell counter vial PS, dust-free capped, 25ml, Ø32x60 mm, suitable for Coulter counter and Hycel counter, tube transparent, cap natural
<b>Dimensions</b>	diameter: 32,00 mm (+/- 0,2 mm)
	height: 60,00 mm (+/- 0,2 mm)
<b>Colour</b>	tube transparent, cap natural
<b>Nominal Capacity</b>	25ml
<b>Material Characteristics</b>	polystyrene (PS)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
<b>Packaging</b>	pcs per bag: 250
	pcs per carton: 1.000
	pcs per pallet: 20.000
<b>Minimum order qty.</b>	1.000 pcs
<b>Storage</b>	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
<b>Approved by nerbe plus quality department</b>	