

HUBERLAB. AG
Industriestrasse 123
4147 Aesch

T 061 717 99 77
F 061 711 93 42

www.huberlab.ch
info@huberlab.ch

HUBERLAB.

committed to science



HETTICH BENELUX
ASPIRATION SYSTEMS



ENGINEERED
FOR YOU

Quality, Innovation and Service

Hettich Benelux is your leading supplier of laboratory equipment, which is tailored to your specific applications and requirements. Expect quality, innovation and service from Hettich Benelux, in everything we do. You will see it in our products, A-class parts, certifications, eco-awareness, service and communication. Not only is innovation important, the maintenance of your equipment is also in good hands.

ASPIRATION SYSTEMS



High ease of use

Quick-lock-adaptors serially for faster bottle removal.

Small device footprint.

Space saving

Aspirators for all applications with Luer-Lock adapters.

Application diversity

Individual flow control and level control.

Secure aspiration

Multifunctional Aspiration Systems AA 02 / 04



Components

Automatic membrane vacuum pump (17 l / min air, 100 mbar abs.)
2 l or 4 l bottle incl. Quick-lock-adaptors, filter, ventilation insert &
3 ml silicon hose
Adjustable handle
Dimensions: 170 x 220 x 155 mm
Electric power: 100 W
Power supply: 115/230V, 50/60 Hz

Article no.	Article description
-------------	---------------------

981323101	AA 02 Aspiration system with 2 l bottle
981323001	AA 04 Aspiration system with 4 l bottle

- High aspiration speed and final vacuum
- Standard Luer-Lock connections

Multifunctional Aspiration Systems A 24



Standard for 2 users: components AA 24

Automatic membrane vacuum pump (17 l / min air, 100 mbar abs.)
4 l bottle incl. Quick-lock-adaptors, filter, ventilation insert &
3 ml silicon hose
2 adjustable handles
Y-distributor
Dimensions: 170 x 220 x 155 mm
Electric power: 100 W
Power supply: 115/230V, 50/60 Hz

Article no.	Article description
-------------	---------------------

981323201	AA 24 Aspiration system with 4 l bottle, 2 handles
-----------	--

- Automatic vacuum pump
- Two users one price

Standard Aspiration Systems AC 02 / 04



25 l/min
(Air)

250 mbar
abs.

Components

Membrane vacuum pump (25 l / min air, 250 mbar abs.)
2 l or 4 l bottle incl. Quick-lock-adapters, filter, ventilation insert &
3 ml silicon hose
Dimensions: 170 x 220 x 155 mm
Electric power: 100 W
Power supply: 115/230V, 50/60 Hz

Article no. Article description

981423101	AC 02 Aspiration system with 2 l bottle
981423001	AC 04 Aspiration system with 4 l bottle

- **Comfortable: Quick-lock-adaptors**
- **Space saving: Small device footprint**

Compact Aspiration Systems AZ 02 / 04



8 l/min
(Air)

700 mbar
abs.

Components

Membran vacuum pump (8 l / min air, 700 mbar abs.)
2 l or 4 l bottle incl. Quick-lock-adapters, filter, ventilation insert &
3 ml silicon hose
Dimensions Pump: 75 x 170 x 60 mm
Electric power: 5 W
Power supply: 115/230V, 50/60 Hz

Article no. Article description

981523101	AZ 02 Aspiration system with 2 l bottle
981523001	AZ 04 Aspiration system with 4 l bottle

- **Ease of use: Quick-lock-adaptors**
- **Very quiet vacuum pump**

Aspiration Devices AE & AK for vacuum connection in the laboratory



Components simple aspiration device

- 2 l or 4 l bottle incl. Quick-lock-adapters, filter, ventilation insert & 3 ml silicon hose


Components comfortable aspiration device

- 2 l or 4 l bottle incl. Quick-lock-adapters, filter, ventilation insert & 3 ml silicon hose
- Bottle holder
- Handle

Article no.	Article description
800015400	AE 02 Aspiration system with 2 l bottle
800015500	AE 04 Aspiration system with 4 l bottle
800014800	AK 02 Aspiration system with 2 l bottle
800014900	AK 04 Aspiration system with 4 l bottle

- **Enjoy the benefits without the pump**
- **Just connect the unit to an existing vacuum system**

Extend the multifunctional aspiration systems up to four workplaces and reduce the costs per user.



Handle and aspiration adapter




90° adapter




Article no.
100210057

8-channel-stainless adapter



Article no.
100210059

4-channel-stainless adapter



Article no.
100210060

More Adapters

Article no.	Article description
100210058	Plastic adapter for Pasteur pipettes and 8-channel
800013700	Adapter with aspiration rake, 8-channel, autoclavable
800015600	Aspiration rake, 8-channel, as spare part
800013300	Stainless steel suction adapter, length 60 mm
800013500	Stainless steel suction adapter, length 200 mm
800013400	Stainless steel suction adapter, length 150 mm





Article no.	Article description
800014600	Spare bottle 2 l with fast lock clutches
800014700	Spare bottle 4 l with fast lock clutches
800014550	Acoustical level control sensor for 2l aspiration system
800014450	Acoustical level control sensor for 4l aspiration system
800010102	Bottle holder for 2 l bottle
800010202	Bottle holder for 4 l bottle
800014400	Cover for bottle 2 l with fast lock clutches
800014500	Cover for bottle 4 l with fast lock clutches
100210028	Hose clip for fast lock clutches
800015800	Safety filter (packing unit = 3 pieces)
800015900	Y-Distributor for connecting several collection unit
800050700	Stainless steel holder for handle
800014000	Foot switch for switching the pump on and off
800015300	AP 04 Membrane vacuum pump, 700mbar

Article no.	Article description
980000101	Country Set Europe
980000201	Country Set Switzerland
980000301	Country Set UK

Vacuum pumps

AP 15 Membrane vacuum pump, 250 mbar.
 Bottle holder is optional.



Article no.
800015000

AP 30 Membrane vacuum pump, 100 mbar.
 Automatic
 Bottle holder is optional.



Article no.
800015100