

Nitrile gloves PRO.TECT



Disposable nitrile glove orange, powder free, extra strong and with Diamondgrip technology

Properties

- Made of synthetic nitrile rubber
- Without silicone
- PRO.TECT Orange HD gloves are thicker and heavier than standard examination gloves
- Made with a full palm and raised fingers Grip pattern on both sides
- Tested according to EN374
- AQL 1.5
- Latex free
- Without powder



Benefits

- Provides improved resistance against snagging and abrasion with maximum dexterity
- Heavy-duty construction is rugged, durable
- Grip pattern provides secure grip, especially in wet or oily conditions
- Suitable for people with latex allergies (Type 1 sensitization)
- Minimizes particle contamination



Order information

Art.-Nr.	Product description	Piece / carton
3.9900.01	Size S 5 - 6, orange	100
3.9900.02	Size M 7 - 8, orange	100
3.9900.03	Size L 8 - 9, orange	100
3.9900.04	Size XL 9 - 10, orange	100

Marked items are included in the sample pack.