

## Product Specificatio

# BLUE PEARL

## Disposable Nitrile Examination Gloves

- Powder free
- Blue
- Chemical resistance protective gloves
- Latex-free
- Suitable for handling food
- Non-sterile



SAFETY MEETS QUALITY

CE 0197



### Material

100% Synthetic Nitrile Latex

### Performance requirements

Characteristics	Inspection Level	Acceptable Quality Level
Dimensions	S-2	4.0
Physical Properties	S-2	4.0
Freedom from Holes	G-1	1.5
Major Visual defects	G-1	1.5
Minor Visual defects	G-1	4.0

### Dimensions

Number	8801	8802	8803	8804	8805
Size	XS = 5-6	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	75 ± 5 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	≥ 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

### Thickness

Measuring Points	Single Wall (mm)
Cuff (distance 25 mm to the beaded edge)	0.08 ± 0.03
Palm (in the middle of the palm)	0.10 ± 0.03
Fingertip (distance 15 mm to top of middle finger)	0.12 ± 0.03

### Physical Properties

Standard	Before Aging	After Aging
Tensile Strength (N)	min. 6.0	min. 6.0

### Powder

No powder lubricant added

### Protein content

No proteins contained

### Cuff

Beaded cuff

### Design

Ambidextrous

### Surface

Micro-roughened finger tips

### Color

Blue

### Quality assurance

G-construction sample test TÜV Rheinland – CE 0197

PSA – EN 374 CAT III, EN 388, EN 420, 89/686/EEC

MPG – EN 455 Teil 1-4

EG-Standard 596/2009

EG-Standard 202/2014

EN ISO 9001:2008

EN ISO 13485:2012

### Packing Style

100 pcs per dispenser, 10 dispensers per carton

### Hint

The above information apply to the standard version.  
Subject to change.