

DURAN® MULTIPORT CONNECTOR CAPS GL 25

committed to
brands.



HUBERLAB.

committed to science



**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN® MULTIPOINT CONNECTOR CAPS GL 25

2, 3 or 4 port connector caps for GL 25 DURAN® bottles


Applications

- Low volume sampling
- Media supply for small bioreactors
- Fed-batch bioreactor
- Perfusion circuit for cell culture
- Oxygenation / aeration of very small samples
- Use in a chemistry micro-reactor

Features

- Caps have a GL 25 thread that fit the smaller 10 or 25 mL DURAN® bottles
- Barbed male hose connector are suitable for tubing with a 3.2 mm (1/8") inner diameter (such as PharMed® BPT Tubing*)
- Available with 2, 3, or 4 surgical grade stainless steel (316L) barbed tubing ports
- Construction materials are certified as approved for food contact
- High temperature resistance (up to 180 °C)
- Autoclavable (@ 121 °C / 15 minutes)
- Reusable with a long service life

*PharMed® is a registered trademark of Saint-Gobain



Ordering information

	Description	Pack Unit	Cat. No.
NEW	DURAN® 2-port Connector caps GL 25, red PBT screw cap, PTFE insert and 2 stainless steel barbs	1	29 260 13 01
NEW	DURAN® 3-port Connector caps GL 25, red PBT screw cap, PTFE insert and 3 stainless steel barbs	1	29 261 13 02
NEW	DURAN® 4-port Connector caps GL 25, red PBT screw cap, PTFE insert and 4 stainless steel barbs	1	29 262 13 03

Recommended compatible DURAN® bottles

Description	Capacity (mL)	Pack Unit	Cat. No.
DURAN® Clear GL 25 bottle supplied without screw cap. Retrace coded.	10	10	21 801 08 02
DURAN® Amber GL 25 bottle supplied without screw cap. Retrace coded.	10	10	21 806 08 07
DURAN® Clear protect safety coated GL 25 bottle supplied without screw cap. Retrace coded.	10	10	21 805 08 06
DURAN® Clear GL 25 bottle supplied without screw cap. Retrace coded.	25	10	21 801 14 04
DURAN® Amber GL 25 bottle supplied without screw cap. Retrace coded.	25	10	21 806 14 09
DURAN® Clear protect safety coated GL 25 bottle supplied without screw cap. Retrace coded.	25	10	10 926 76
DURAN® Amber protect safety coated GL 25 bottle supplied without screw cap. Retrace coded.	25	10	21 806 14 33

These products are available from your specialist laboratory dealer.



Contact.

HUBERLAB.

committed to science

OneStop. OneShop – HUBERLAB. bietet das Rundum-Sorglos-Paket für alle Bedürfnisse im Labor.

Profitieren Sie von unseren langjährigen Partnerschaften zu Top-Lieferanten.

Als Vollversorger für das Labor lässt unser umfassendes Sortiment keine Wünsche offen.

OneStop. OneShop – HUBERLAB. Nous offre le package sans souci pour tous les besoins de votre laboratoire.

Profitez de nos partenariats de longue date avec les meilleurs fournisseurs.

En tant que fournisseur complet pour le laboratoire, notre gamme complète ne laisse rien à désirer.

OneStop. OneShop – HUBERLAB. offers the all-included service package for all of your laboratory requirements.

Benefit from our long-standing partnerships with top-suppliers.

As a full-range supplier for the laboratory, our comprehensive range leaves nothing to be desired.



Wir beraten Sie gerne persönlich
Nous vous conseillons personnellement
We'll be happy to advise you personally



T +41 61 717 99 77



F +41 61 711 93 42



info.huberlab.ch



www@huberlab.ch



Let's chat
www.huberlab.ch



HUBERLAB. AG
Industriestrasse 123
4147 Aesch