

# DURAN® MULTIPOINT CONNECTOR CAPS GL 32

committed to  
**brands.**



**HUBERLAB.**

committed to science



**DURAN  
WHEATON  
KIMBLE**

Excellence in your hands

# DURAN® MULTIPOINT CONNECTOR CAPS GL 32

2, 3 or 4 port connector caps for GL 32 DURAN® bottles

## Applications

- Low volume sampling
- Media supply for small bioreactors
- Fed-batch bioreactor
- Perfusion circuit for cell culture
- Oxygenation / aeration of very small samples
- Use in a chemistry micro-reactor

## Features

- Caps have a GL 32 thread that fit the smaller 50 mL round or 100 mL DURAN® square bottles
- Barbed male hose connector are suitable for tubing with a 3.2 mm (1/8") inner diameter (such as PharMed® BPT Tubing\*)
- Available with 2, 3, or 4 surgical grade stainless steel (316L) barbed tubing ports
- Construction materials are certified as approved for food contact
- High temperature resistance (up to 180 °C)
- Autoclavable (@ 121 °C / 15 minutes)
- Reusable with a long service life

\*PharMed® is a registered trademark of Saint-Gobain



## Ordering information

	Description	Pack Unit	Cat. No.
NEW	DURAN® 2-port Connector caps GL 32, red PBT screw cap, PTFE insert and 2 stainless steel barbs	1	29 260 19 01
NEW	DURAN® 3-port Connector caps GL 32, red PBT screw cap, PTFE insert and 3 stainless steel barbs	1	29 261 19 02
NEW	DURAN® 4-port Connector caps GL 32, red PBT screw cap, PTFE insert and 4 stainless steel barbs	1	29 262 19 03

## Recommended compatible DURAN® bottles

Description	Capacity (mL)	Pack Unit	Cat. No.
DURAN® Laboratory Bottle GL 32, clear, without screw cap. Retrace coded.	50	10	21 801 17 04
DURAN® Laboratory Bottle GL 32, amber, without screw cap. Retrace coded.	50	10	21 806 17 09
DURAN® Protect Laboratory Bottle GL 32, clear, safety coated, without screw cap. Retrace coded.	50	10	10 926 77
DURAN® Laboratory Bottle Square GL 32, clear, without screw cap. Retrace coded.	100	10	21 820 24 04

These products are available from your specialist laboratory dealer.



# Contact.

## HUBERLAB.

committed to science

**OneStop. OneShop – HUBERLAB. bietet das Rundum-Sorglos-Paket für alle Bedürfnisse im Labor.**

Profitieren Sie von unseren langjährigen Partnerschaften zu Top-Lieferanten.

Als Vollversorger für das Labor lässt unser umfassendes Sortiment keine Wünsche offen.

---

**OneStop. OneShop – HUBERLAB. Nous offre le package sans souci pour tous les besoins de votre laboratoire.**

Profitez de nos partenariats de longue date avec les meilleurs fournisseurs.

En tant que fournisseur complet pour le laboratoire, notre gamme complète ne laisse rien à désirer.

---

**OneStop. OneShop – HUBERLAB. offers the all-included service package for all of your laboratory requirements.**

Benefit from our long-standing partnerships with top-suppliers.

As a full-range supplier for the laboratory, our comprehensive range leaves nothing to be desired.



**Wir beraten Sie gerne persönlich**  
Nous vous conseillons personnellement  
We'll be happy to advise you personally



T +41 61 717 99 77



F +41 61 711 93 42



info.huberlab.ch



www@huberlab.ch



Let's chat  
www.huberlab.ch



**HUBERLAB. AG**  
Industriestrasse 123  
4147 Aesch