

committed to
brands.



**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN® PURE

DEVELOPED FOR PHARMA

DURAN® PURE is a range of high-quality bottles and closures developed for the needs of the pharmaceutical industry. **DURAN® PURE** meets the requirements of the international pharmacopeias and is available with the certificates required by the regulatory authorities.

FEATURES

- **EP/USP Type I borosilicate glass**
- **Size range 25 ml - 20 l**
- **Protective cover** prevents contamination
- **Single batch delivery** CoA available
- **100% traceability**
- **100% visual inspection**
- **Change control**
- **Drug Master File**

APPLICATION AREAS

- Manufacturing of pharmaceutical products
- Storage of bulk drug intermediates
- Transportation of APIs



DURAN® PURE Bottles

Article number clear glass	Article number amber coloured	Volume [ml]	DIN thread GL	Bottles per box [pcs.]	Bottles per pallet [pcs.]
81 801 14 01	81 806 14 06	25	25	10	940
81 801 24 06	81 806 24 02	100	45	10	540
81 801 36 08	81 806 36 04	250	45	10	280
81 801 44 07	81 806 44 03	500	45	10	360
81 801 54 03	81 806 54 08	1 000	45	10	210
81 801 63 05	81 806 63 01	2 000	45	10	120
81 801 73 01	81 806 73 06	5 000	45	1	48
81 801 86 06	81 806 86 02	10 000	45	1	30
81 801 91 05	81 806 91 01	20 000	45	1	12

ALSO AVAILABLE: DURAN® PURE GLS 80® bottles and screw caps

DURAN® PURE Premium Screw Caps

DURAN® PURE Premium Screw Caps made from PFA (fluoropolymer) are chemically and thermally very resistant and compatible with autoclaving and other common sterilization methods. The products include silicone seals which are coated with highly inert PTFE.

FEATURES

- Caps made from PFA with PTFE coated silicone seal
- Produced and packaged in clean rooms
- Single batch delivery CoA available
- 100% traceability
- Change control
- Drug Master File



DURAN® PURE Premium Screw Cap



DURAN® PURE Premium Lip Seal Cap

Article number	DIN thread GL	Description	Screw Caps per PE bag [pcs.]
89 108 13 06	25	DURAN® PURE Premium Screw Cap	5
89 108 86 79	45	DURAN® PURE Premium Screw Cap	5
89 886 28 02	45	DURAN® PURE Premium Lip Seal Cap	5



HUBERLAB.
committed to science

Contact.

HUBERLAB.

committed to science

OneStop. OneShop – HUBERLAB. bietet das Rundum-Sorglos-Paket für alle Bedürfnisse im Labor.

Profitieren Sie von unseren langjährigen Partnerschaften zu Top-Lieferanten.

Als Vollversorger für das Labor lässt unser umfassendes Sortiment keine Wünsche offen.

OneStop. OneShop – HUBERLAB. Nous offre le package sans souci pour tous les besoins de votre laboratoire.

Profitez de nos partenariats de longue date avec les meilleurs fournisseurs.

En tant que fournisseur complet pour le laboratoire, notre gamme complète ne laisse rien à désirer.

OneStop. OneShop – HUBERLAB. offers the all-included service package for all of your laboratory requirements.

Benefit from our long-standing partnerships with top-suppliers.

As a full-range supplier for the laboratory, our comprehensive range leaves nothing to be desired.



Wir beraten Sie gerne persönlich
Nous vous conseillons personnellement
We'll be happy to advise you personally



T +41 61 717 99 77



F +41 61 711 93 42



info.huberlab.ch



www@huberlab.ch



Let's chat
www.huberlab.ch



HUBERLAB. AG
Industriestrasse 123
4147 Aesch