

committed to
brands.



**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN® PURE GLS 80®

DEVELOPED FOR PHARMA

DURAN® PURE is a range of high-quality bottles and closures developed for the needs of the pharmaceutical industry. DURAN® PURE meets the requirements of the international pharmacopeias and is available with the certificates required by the regulatory authorities.

FEATURES

- EP/USP Type I borosilicate glass
- Size range 500 ml - 20 l
- GLS 80® thread for fast closing & opening
- Protective cover prevents contamination
- Single batch delivery CoA available
- 100% traceability
- 100% visual inspection
- Change control
- Drug Master File

APPLICATION AREAS

- Manufacturing of pharmaceutical products
- Storage of bulk drug intermediates
- Transportation of APIs



DURAN® PURE GLS 80® Bottles

Article number clear glass	Article number amber coloured	Volume [ml]	Bottles per box [pcs.]	Bottles per pallet [pcs.]
81 860 44 03	81 866 44 09	500	10	280
81 860 54 08	81 866 54 05	1 000	10	210
81 860 63 01	81 866 63 07	2 000	10	120
81 860 73 06	81 866 73 03	5 000	1	48
81 860 86 02	81 866 86 08	10 000	1	30
81 860 91 01	81 866 91 07	20 000	1	12

ALSO AVAILABLE: DURAN® PURE bottles and screw caps with GL 45 thread

DURAN® PURE GLS 80® Screw Cap

The DURAN® PURE GLS 80® Screw Cap made from PSU is chemically and thermally very resistant and compatible with autoclaving and other common sterilization methods. The product includes a silicone seal which is coated with highly inert PTFE.

FEATURES

- Cap made from PSU with PTFE coated silicone seal
- GLS 80® thread for fast closing & opening
- Single batch delivery CoA available
- 100% traceability
- Change control



Article number	Screw Caps per PE bag [pcs.]	Screw Caps per pallet [pcs.]
89 116 91 09	1	1435



HUBERLAB.
committed to science

Contact.

HUBERLAB.

committed to science

OneStop. OneShop – HUBERLAB. bietet das Rundum-Sorglos-Paket für alle Bedürfnisse im Labor.

Profitieren Sie von unseren langjährigen Partnerschaften zu Top-Lieferanten.

Als Vollversorger für das Labor lässt unser umfassendes Sortiment keine Wünsche offen.

OneStop. OneShop – HUBERLAB. Nous offre le package sans souci pour tous les besoins de votre laboratoire.

Profitez de nos partenariats de longue date avec les meilleurs fournisseurs.

En tant que fournisseur complet pour le laboratoire, notre gamme complète ne laisse rien à désirer.

OneStop. OneShop – HUBERLAB. offers the all-included service package for all of your laboratory requirements.

Benefit from our long-standing partnerships with top-suppliers.

As a full-range supplier for the laboratory, our comprehensive range leaves nothing to be desired.



Wir beraten Sie gerne persönlich
Nous vous conseillons personnellement
We'll be happy to advise you personally



T +41 61 717 99 77



F +41 61 711 93 42



info.huberlab.ch



www@huberlab.ch



Let's chat
www.huberlab.ch



HUBERLAB. AG
Industriestrasse 123
4147 Aesch