

FACEMATE (18.2021.02)

Flat Fold Respirator Masks

Lightweight design • white • Breathe Easy in Comfort •
Adjustable nose piece • latex free •
Filter efficiency 94%

Colour: white

Item Number and Colour:

18.2021.02 | White

Packaging:

25 pieces per Box
4 Boxes per Carton

Classification:



EN 149:2001 + A1:2009



Product Info

To protect people working in environments where harmful particles, pathogens and aerosols present a danger to their respiratory systems and long-term health. Healthcare workers, Medical and Military personnel, Cleanroom operatives, Industrial, Agricultural and Construction workers will benefit from using the new Facemate Respirator.

Dimensions

Item-Nr	Filter efficiency	Breathin resistance (mbar)30L/ min	Defined protection Level	Nominal protection Factor
18.2021.02	94% to 0.6µm	< 0.70	10 OEL	12 OEL

Comfort Protection

- Horizontal fold style for greater Protection and Comfort
- Individually wrapped
- Tested for a wide range of face profiles and head sizes
- Meets FFP2 Classification in accordance with EN 149:2001 + A1: 2009
- Latest filter design technology providing lowest possible breathing resistance
- Hypoallergenic materials providing extra comfort during wearing for extended periods
- Maximum wearing time of 8 hours

Hint

The above information applies to the standard version.
Subject to changes

