



Your **Power** for Health

greiner bio-one

Order Form

Barcode Service for Cryo.s™ Datamatrix Racks and Storage Boxes

Customer No.		Phone	
Name		Fax	
Institution		E-Mail	
Street		PO # / Date	
ZIP / City		Shipping Address (if different from the left)	
Country			

① Please choose **ONE** of the following products by ticking one box below!
(For ordering further products please fill in a new Order Form).

Order no.	803202-128	803270-128	802576-128	802506-128
Capacity	48 tubes	48 tubes	81 tubes	81 tubes
Storable tube sizes	1 / 2 ml	4 ml	1 / 2 ml	4 ml
Rack height [mm]	52.5	88.5	52.0	88.1
Quantity per bag / case	20 / 20	15 / 15	1 / 10	1 / 10
Tick here	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Datamatrix Code	<p>ECC200 14x14 dots on rack bottom below wells A1, A2, B1 and B2 Bottom view</p>	<p>ECC200 14x14 dots on rack bottom below wells Bottom view</p>
	<p>Barcode type 128 on S1 side of rack Side view</p>	<p>Barcode type 128 on rack side Front view</p>
Linear Barcode		

② Please fill in the order quantity and code sequences!

Order quantity in pcs.	Code sequences
	First code:
	Last code:

Please note: Code length must not exceed the maximum data capacity of 16 digits for pure numeric codes or 10 characters for codes composed of alphanumeric characters.

<p>Order released by:</p> <p>Company Name: _____</p> <p>Customer Name: _____</p> <p>Date / Signature: _____</p>
--