

LaboPlast® Laboratory stands

NEW

- ▶ Funnel stands /
Filtering stands
- ▶ Stand for Imhoff
sedimentation funnel
- ▶ Pipette stand

With innovative
magnet concept

- ▶ Easy and flexible
height adjustment
- ▶ Comfortable and
smooth operation



LaboPlast

Funnel stands

Sturdy and stable funnel stand, resistant to chemicals, shatter-proof and with continuously variable height adjustment.

The funnel stand also makes an excellent filtering frame.

Thanks to an innovative magnet concept no clamping screws are necessary. The height of the funnel holder can be adjusted easily and flexibly with one hand.

The funnel stand is suitable for all funnels from 40 to 180 mm in diameter. The funnel holders can be used with reducing inserts for small funnels.

- ▶ PP, chrome-plated steel
- ▶ Continuously variable height adjustment with innovative magnet concept
- ▶ Adjustable height up to 450 mm
- ▶ For funnels Ø 40 - 180 mm
- ▶ Additional funnel holders are available as accessories



Funnel stand for two funnels

Funnel stand with two funnel holders and two reducing inserts. Space-saving design with two practical storage areas. The lab frame can be flexibly adapted with corresponding accessories, for example with adjustment funnel holders or a holder for Imhoff funnels.



Funnel stand for two funnels

Model	Item no.
for two funnels	9622-0200 N4
Accessories	
Description	Item no.
① Funnel holder with reducing insert, for funnels Ø 40-180 mm	9622-0001 N4
② Sedimentation funnel holder for Imhoff funnels from Ø 90 mm	7118-2001 N4

Funnel stand for four funnels

Funnel stand with four funnel holders and four reducing inserts. A very sturdy pedestal provides a storage area for four containers and ensures absolute safety against overturning. The lab frame can be flexibly adapted with corresponding accessories, for example with adjustment funnel holders or a holder for Imhoff funnels. The funnel stand can also be converted into a pipette stand.



Funnel stand for four funnels

Model	Item no.
for four funnels	9622-0400 N4
Accessories	
Description	Item no.
① Funnel holder with reducing insert, for funnels Ø 40-180 mm	9622-0001 N4
② Sedimentation funnel holder for Imhoff funnels from Ø 90 mm	7118-2001 N4
③ Replacement set for pipette stand, two perforated plates/one drip pan	9630-0100 N4

Stand for Imhoff sedimentation funnel

The stand for sedimentation funnels is used to keep, hold or store Imhoff funnels safely and reliably during sedimentation.

Two funnel holders ensure that Imhoff funnels will be held safely and securely and can be placed in an absolutely vertical position.

Stand with two holders for sedimentation funnels, suitable for funnels 90 mm or more in diameter.

Thanks to an innovative magnet concept no clamping screws are necessary. The height of the funnel holder can be adjusted easily and flexibly with one hand.

With the corresponding accessory the holders for Imhoff funnels can also be replaced by holders for all-purpose funnels.

- ▶ Stand for Imhoff funnel
- ▶ PP, chrome-plated steel
- ▶ Adjustable height up to 450 mm
- ▶ For funnels with Ø 90 mm or more



Stand for Imhoff sedimentation funnel

Model	Item no.
for two funnels	7118-2000 N4

Accessories

Description	Item no.
① Funnel holder with reducing insert, for funnels Ø 40-180 mm	9622-0001 N4
② Sedimentation funnel holder for Imhoff funnels from Ø 90 mm	7118-2001 N4

Pipette stand



Frame for pipettes with two rotating perforated plates and drip pan. The holders are suitable for up to 40 pipettes of different sizes, for all measuring and full pipettes up to a max. diameter of 18 mm.

A drip pan protects the tips of the pipettes and collects liquid as it drips off.

The pipette stand is suitable for drying and for space-saving storage of pipettes or thermometers, test tubes and similar items. A very sturdy pedestal ensures absolute safety against overturning.

No clamping screws are necessary thanks to an innovative magnet concept. The height of the holding plates can be adjusted easily and flexibly with one hand.

The lab frame can be flexibly adapted with corresponding accessories, for example as a frame for Imhoff sedimentation funnels or as a funnel frame.

- ▶ PP, chrome-plated steel
- ▶ Adjustable height up to 450 mm
- ▶ 40 slots, up to Ø 18 mm



Pipette stand

Model	Item no.
for 40 pipettes	9630-0090 N4

Accessories

Description	Item no.
① Funnel holder with reducing insert, for funnels Ø 40-180 mm	9622-0001 N4
② Sedimentation funnel holder for Imhoff funnels from Ø 90 mm	7118-2001 N4