

Technical Datasheet

Date: 13.10.17

Page 1 of 1

		EN ISO certified
Article - No:	02-322-1056	Revision: 1
Remarks		single use only

Specification

Product	Test tube PP with enclosed screw cap PE, 5ml, Ø16x54 mm, flat bottom, tube transparent, cap white, max. RCF 3.000g, Autocl. up to 121°C (no cap)
Dimensions	diameter: 16,00 mm (+/- 0,2 mm)
	height (without screw cap): 54,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap white
Nominal Capacity	5ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	cap: -80°C up to +75°C
	container: -20°C up to +121°C
Autoclavability	Autocl. up to 121°C (no cap)
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 500
	pcs per carton: 2.500
	pcs per pallet: 50.000
Storage	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	