

Technical Datasheet

Date: 15.10.18

Page 1 of 1

		EN ISO certified
Article - No:	11-411-6000	Revision: 2
Remarks		single use only

Specification

Product	Screw cap container PP, 20ml, container transparent, cap green, Moulded rings every 10 ml, grading on 3 sites, sterile EO
Colour	container transparent, cap green
Graduation	Moulded rings every 10 ml, grading on 3 sites
Nominal Capacity	20ml
Material Characteristics	container: polypropylene (PP) cap: high-density polyethylene (HDPE)
Autoclavability	autocl. up to 121°C
Purity Grade	do not contain any pyrogens
Sterilisation Process	gas sterilisation by ethyleneoxid
Shelf - Life	24 months at controlled storage conditions
Packaging	pcs per carton: 600 pcs per pallet: 18.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	