

# Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>11-622-0000</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Screw cap container PP, 20ml, Ø31x48 mm, with enclosed screw cap PE, container transparent, cap natural
<b>Dimensions</b>	diameter: 31,00 mm (+/- 0,2 mm)
	height: 48,00 mm (+/- 0,2 mm)
<b>Colour</b>	container transparent, cap natural
<b>Nominal Capacity</b>	20ml
<b>Material Characteristics</b>	container: polypropylene (PP)
	cap: polyethylene (PE)
<b>Temperature range</b>	-10°C - 40°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	flexible and unbreakable
<b>Shelf - Life</b>	not limited at controlled storage conditions
<b>Packaging</b>	pcs per bag: 1.000
	pcs per carton: 1.000
	pcs per pallet: 24.000
<b>Minimum order qty.</b>	1.000 pcs
<b>Storage</b>	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
<b>Approved by nerbe plus quality department</b>	