

Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-652-0000	Revision: 1
Remarks		single use only

Specification

Product	Screw cap container PP, 100ml, Ø48x80 mm, with enclosed screw cap PE, container transparent, cap natural
Dimensions	diameter: 48,00 mm (+/- 0,2 mm)
	height: 80,00 mm (+/- 0,2 mm)
Colour	container transparent, cap natural
Nominal Capacity	100ml
Material Characteristics	container: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-10°C - 40°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	flexible and unbreakable
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 250
	pcs per carton: 250
	pcs per pallet: 6.000
Minimum order qty.	500 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	