

Technical Datasheet

Date: 17.11.16

Page 1 of 1

EN ISO certified		
Article - No:	11-712-1010	Revision: 1
Remarks	suitable for formalin and chloroform	single use only

Specification

Product	Screw cap container PP, 1.250ml, Ø106x159 mm, with fitted screw cap PE, with EPE-pad, white, autoclavable up to 121°C (without pad)
Colour	white
Nominal Capacity	1.250ml
Material Characteristics	container: polypropylene (PP) cap: polypropylene (PP) pad: expanded polyethylene (EPE)
Temperature range	-80°C - 70°C
Autoclavability	autoclavable up to 121°C (without pad)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004 flexible and soft
Packaging	pcs per cs: 35 pcs per carton: 35 pcs per pallet: 700
Minimum order qty.	35 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	