

Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-802-1150	Revision: 1
Remarks		single use only

Specification

Product	Multi purpose container PP, 200ml, Ø66x88 mm, with enclosed snap-on lid PE, transparent, graduated
Dimensions	outer diameter: 66,00 mm (+/- 0,2 mm)
	inner diameter: 62,00 mm (+/- 0,2 mm)
	height (with lid): 88,00 mm (+/- 0,2 mm)
	height (without lid): 87,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	graduated
Nominal Capacity	200ml
Material Characteristics	container: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-80°C - 75°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Packaging	pcs per bag: 660
	pcs per carton: 660
	pcs per pallet: 10.560
Minimum order qty.	660 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	