

Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-842-6100	Revision: 1
Remarks		single use only

Specification

Product	Multi purpose container PP, 1.000ml, Ø110x132 mm, with fitted snap-on lid PE , transparent, graduated, sterile R
Dimensions	outer diameter: 110,00 mm (+/- 0,2 mm)
	inner diameter: 107,00 mm (+/- 0,2 mm)
	height (with lid): 132,00 mm (+/- 0,2 mm)
	height (without lid): 131,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	graduated
Nominal Capacity	1.000ml
Material Characteristics	container: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-80°C - 75°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 80
	pcs per carton: 80
	pcs per pallet: 960
Minimum order qty.	80 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	