

Technical Datasheet

Date: 17.11.16

Page 1 of 1

		EN ISO certified
Article - No:	11-852-0100	Revision: 1
Remarks		single use only

Specification

Product	Multi purpose container PP, 30ml, Ø38x42 mm, transparent, graduated, autoclavable up to 121°C
Dimensions	diameter: 38,00 mm (+/- 0,2 mm)
	height: 42,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	graduated
Nominal Capacity	30ml
Material Characteristics	polypropylene (PP)
Temperature range	-196°C - 121°C
Autoclavability	autoclavable up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 75
	pcs per carton: 3.750
	pcs per pallet: 180.000
Minimum order qty.	3.750 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	