

Technical Datasheet

Date: 14.02.17

Page 1 of 1



EN ISO 13485 certified

		EN ISO certified
Article - No:	02-162-0001	Revision: 8
Remarks	Fitting plug: 02-633-0000	single use only

Specification

Product	Test tube PP, round bottom, 5ml, Ø13x75 mm, transparent, max. RCF 3.000g, autoclavable up to 121°C
Dimensions	outer diameter: 13,00 mm (+/- 0,2 mm)
	height: 75,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	5ml
Material Characteristics	polypropylene (PP)
Temperature range	-10°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 4.000
	pcs per pallet: 80.000
Minimum order qty.	4.000 pcs
Storage	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	