

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	02-352-1000	Revision: 1
Remarks		single use only

Specification

Product	Test tube PP with enclosed screw cap PE, 12ml, Ø16x100 mm, round bottom, tube transparent, cap natural, graduated, sterile, max. RCF 4.000g
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm)
	height (withour screw cap): 100,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Graduation	graduated
Nominal Capacity	12ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-20°C - 65°C
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert
	physiologically harmless
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 1.000
	pcs per pallet: 24.000
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	