

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	02-352-6000	Revision: 2
Remarks		single use only

Specification

Product	Test tube PP with fitted screw cap PE , 12ml, Ø16x100 mm, round bottom, tube transparent, cap natural, graduated, sterile R, max. RCF 4.000g
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm)
	height (withour screw cap): 100,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Graduation	graduated
Nominal Capacity	12ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	tube: -196°C up to +121°C
	cap: -50°C up to +80°C
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 1.000
	pcs per pallet: 24.000
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	