

# Technical Datasheet

Date: 13.12.17

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>02-801-2000</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Cell counter vial PS, dust-free capped, 15ml, Ø22x45 mm, suitable for Labora, Contraves Hellige, Brand, tube transparent, cap natural
<b>Dimensions</b>	diameter: 22,00 mm (+/- 0,2 mm)
	height: 45,00 mm (+/- 0,2 mm)
<b>Colour</b>	tube transparent, cap natural
<b>Nominal Capacity</b>	15ml
<b>Material Characteristics</b>	polystyrene (PS)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
<b>Packaging</b>	pcs per bag: 500
	pcs per carton: 2.000
	pcs per pallet: 24.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	