

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	02-381-8107	Revision: 1
Remarks		single use only

Specification

Product	Test tube PS with fitted screw cap PE, 19ml, Ø16x150 mm, round bottom, tube transparent, cap red, sterile R
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm) height (without screw cap): 150,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap red
Nominal Capacity	19ml
Material Characteristics	tube: polypropylene (PP) cap: polyethylene (PE)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 100 pcs per carton: 800 pcs per pallet: 19.200
Minimum order qty.	800 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	