

# Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>05-733-1000</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Tubes PE for scintillation bottles, 4ml, enclosed screw cap PP, Ø16x57 mm, natural
<b>Dimensions</b>	diameter: 16,00 mm (+/- 0,2 mm)
	height: 57,00 mm (+/- 0,2 mm)
<b>Colour</b>	natural
<b>Nominal Capacity</b>	4ml
<b>Material Characteristics</b>	polyethylene (PE)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 1.000
	pcs per carton: 2.000
	pcs per pallet: 64.000
<b>Minimum order qty.</b>	2.000 pcs
<b>Storage</b>	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
<b>Approved by nerbe plus quality department</b>	