

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	09-703-5000	Revision: 1
Remarks		single use only

Specification

Product	Wire closure bags PE, 540ml, 114x230 mm, transparent, sterile R
Dimensions	width: 114,00 mm (+/- 0,2 mm)
	length: 230,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	540ml
Material Characteristics	Polyvinyl chloride (PVC)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	leak-proof
	odor-proof
Purity Grade	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 500
	pcs per carton: 1.000
	pcs per pallet: 65.000
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	