

# Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
<b>Article - No:</b>	<b>09-723-5000</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Wire closure bags PE, 1.260ml, 140x380 mm, transparent, sterile R
<b>Dimensions</b>	width: 140,00 mm (+/- 0,2 mm)
	length: 380,00 mm (+/- 0,2 mm)
<b>Colour</b>	transparent
<b>Nominal Capacity</b>	1.260ml
<b>Material Characteristics</b>	bag: polyethylene (PE)
	closure: steel wires covered with pressure sensitive tape with a backing of
	polypropylene (PP) or Polyvinyl chloride (PVC)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	leak-proof
	odor-proof
<b>Purity Grade</b>	free of RNase and DNase
	free of pyrogens
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per inner cs: 500
	pcs per carton: 1.000
	pcs per pallet: 65.000
<b>Minimum order qty.</b>	1.000 pcs
<b>Storage</b>	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
<b>Approved by nerbe plus quality department</b>	