

Technical Datasheet

Date: 14.02.17

Page 1 of 1

EN ISO 13485 certified

		EN ISO certified
Article - No:	08-602-0030	Revision: 3
Remarks		single use only

Specification

Product	Universal container PP, 125ml, Ø55x70 mm, with screw thread, without cap, transparent, graduated, with writing area, np pcr ready, autoclavable up to 121°C, CE/IVD
Classification	CE/IVD
Dimensions	upper diameter: 48,00 mm (+/- 0,2 mm) bottom diameter: 55,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	grad., writing area
Nominal Capacity	125ml
Material Characteristics	polypropylene (PP)
Temperature range	-40°C-75°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Packaging	pcs per bag: 750 pcs per carton: 750 pcs per pallet: 18.000
Minimum order qty.	750 pcs
Storage	dry condition: 40%-60% rel. humidity bags & cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	